Detaseal® hydroflow mono

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Detaseal[®] hydroflow mono

Instructions for use

VPS precision impression material, medium bodied consistency – snap set effect

Fig. 1

1. Tray preparation

Select a tray type suited for desired application (standard-size or custom, in metal or hard plastics). To avoid teeth contact with the tray, place stops from wax or light-curing resin into the tray at non-relevant areas of the impression. For optimal adhesion of impression material, we recommend applying a thin layer of **trayloc[®] A** to tray (Fig. 1). Let dry for approx. 3 minutes.

Fig. 2

2. Mixing and dispensing

Insert cartridge into Automix2-System mixing gun (Fig. 2). Remove and dispose of closure cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved (Fig. 3). Attach mixing cannula by aligning notched side to cartridge and turning clockwise until fastened (Fig. 4). Dispenser is now ready to mix silicone in any amount needed. After each use leave the used mixing cannula in place to seal the cartridge.


Fig. 3

3. Procedure

At the beginning introduce a sufficient quantity of **Detaseal[®] hydroflow mono** to the impression tray. Then syringe around the relevant areas with **Detaseal[®] hydroflow lite** or **Detaseal[®] hydroflow Xlite** intraorally. For this purpose intraoral tips or application syringes (Figs. 5 and 6) can be used. Position the loaded tray in the patient's mouth, applying slight pressure. Fix the impression without pressure until the material has set completely (Fig. 7).

Fig. 4

4. Disinfecting

After removing tray from mouth, rinse impression with lukewarm water. If desired, impression can then be disinfected by immersing in 2 % glutaraldehyde for 15 min.

Fig. 5

5. Pouring the model

We recommend waiting at least 30 minutes before pouring the model. The impression will retain its dimensional stability for up to 14 days. All class III and IV dental plasters and standard modelling resins can be used.


Fig. 6

6. Cleaning the tray

After setting, material can be removed with a blunt instrument. Remaining film of **trayloc[®] A** is dissolved by immersion in acetone or isopropanol (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.

Fig. 7


Important working hints

- Substances such as retraction solutions can adversely affect the setting reaction and should therefore be thoroughly rinsed off before taking the impression.
- Do not combine with condensation curing silicones.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- Be sure to remove all traces of impression material from patient's mouth.
- Contact with latex examination gloves or other latex material may affect the setting reaction of **Detaseal[®] hydroflow mono**.
- We recommend use of **Detaseal[®] hydroflow mono lite** and **-Xlite** for double-mix impression techniques and for corrective impressions with two-stage putty and wash techniques.
- In order to enhance easy mouth removal in case of one-step impressions and to avoid blocking in strong undercuts, fill undercuts beforehand with **Detaseal[®] antilock**.
- Please follow the instructions of the safety data sheet!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 072 43/510-0 · Fax: 072 43/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

Made in Germany

08/2017

Indications for use:

- single-mix impression technique
- double-mix impression technique
- situation impressions
- over-ring and fixing impressions
- Impressions for indirect denture relining
- Transfer impressions
- Tray material for the sandwich and double mix impression technique

Technical data:

DIN EN ISO 4823 - Type 2

- Mixing volume:** 50 ml (cartridge)
- Mixing ratio:** 1:1
- Colour code:** base: yellow catalyst: white
- Mixing time:** omitted (Automix2-System)
- Working time:** approx. 1 min. 30 sec.*
- Setting time:** approx. 3 min. 30 sec.*
- Time in mouth:** approx. 2 min.
- Strain in compression:** approx. 3 %
- Recovery from deformation:** approx. 99,8 %
- Linear dimensional change:** approx. 0,2 %
- Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity.
- Storage:**
 - * from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperature accelerate, decreased temperatures retard a. m. times.

Ordering information:

Detaseal[®] hydroflow mono
Standard packing **02754**
2 cartridges of 50 ml
6 mixing cannulas, green

Multi packing **02757**
4 cartridges of 50 ml
12 mixing cannulas, green

Mixing cannulas **02705**
green, 48 pcs.

Detaseal[®] hydroflow lite
Standard packing
2 cartridges of 50 ml
12 mixing cannulas, yellow
regular set **02730**
fast set **02733**

Detaseal[®] hydroflow Xlite
Standard packing
2 cartridges of 50 ml
12 mixing cannulas, yellow
regular set **02741**
fast set **02744**

Detaseal[®] antilock
Standard pack **03535**
1 x 50 ml cartridge
6 Mixing cannulas, yellow
6 Intra-oral tips, yellow

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